

Vivid Vision



HOW IT WORKS

Virtual Reality - This powerful technology allows us to create immersive vision therapy tools that appear to be a life size, 3D virtual environment without the boundaries we are accustomed to with TV and computer screens.



Using a head-mounted display, your patient can be sitting or standing and moving around freely



Work on vergence ranges and jump duction to increase fusion ranges



Full control of elements seen in the amblyopic/strabismic eye



Requires no filters, batteries, or large TV monitors



Start activities at the patient's angle of deviation and decrease the angle as they improve



Distance and near activities all in one device